



ARTHROSCOPY ELECTRODES



- Designed for knee, shoulder, hip, ankle, elbow, & wrist arthroscopy.
- Slim shaft construction allows enhanced access to confined spaces within the joint.
- Each electrode is fitted with an anti rotation collar & can be supplied in a selection of shaft lengths.
- Custom designs available upon request.



ARTHROSCOPY ELECTRODES

Small Joint Arthroscopy Electrodes - 115mm Shaft Wrist - Ankle - Elbow

17-460	Lateral Release	Angled Blade
17-461	Lateral Release	Straight
17-462	Spike	Angled 90°
17-464	Hockey Stick	Angled 90°
17-466	Ball 1 mm	Straight
17-467	Ball 1 mm	Angled 30°
17-468	Ball 2 mm	Straight
17-469	Ball 2 mm	Angled 30°
17-470	Ball 3 mm	Straight
17-471	Ball 3 mm	Angled 30°
17-480	Sickle	Straight

Knee Arthroscopy Electrodes - 150mm Shaft


17-260	Lateral Release	Angled Blade
17-261	Lateral Release	Straight
17-262	Spike	Angled 90°
17-264	Hockey Stick	Angled 90°
17-268	Ball 2 mm	Straight
17-269	Ball 2 mm	Angled 30°
17-270	Ball 3 mm	Straight
17-271	Ball 3 mm	Angled 30°
17-272	Ball 4 mm	Straight
17-273	Ball 4 mm	Angled 30°
17-274	Ball 5 mm	Straight
17-275	Ball 5 mm	Angled 30°
17-280	Sickle	Straight

Shoulder Arthroscopy Electrodes - 190mm Shaft

17-360	Lateral Release	Angled Blade
17-361	Lateral Release	Straight
17-362	Spike	Angled 90°
17-364	Hockey Stick	Angled 90°
17-368	Ball 2 mm	Straight
17-369	Ball 2 mm	Angled 30°
17-370	Ball 3 mm	Straight
17-371	Ball 3 mm	Angled 30°
17-372	Ball 4 mm	Straight
17-373	Ball 4 mm	Angled 30°
17-380	Sickle	Straight

Hip Arthroscopy Electrodes - 250mm Shaft

17-570	Ball 3 mm	Straight
17-571	Ball 3 mm	Angled 30°



Scale 1:1