

# Having trouble fitting your Bipolar Leads into your generators?

Trewavis Surgical has the solution!  
 Introducing a new range of Australian Made ESU Bipolar Adaptors to suit a variety of generator models.

## Bipolar Adaptors




Product Number: 300-770-C  
 To suit: Conmed generators (ABC, 7550, Excalibur)




Product Number: 300-770-D  
 To suit: Erbe generators (ICC80, ICC50) and Alcon Infiniti



Product Number: 300-770-E  
 To suit: Alcon generators (Acurs)



Trewavis Twin Pin Bipolar Lead  
 To suit Valleylab or similar twin pin spacing.



Remember "Flying Leads" are now prohibited for sale in Australia by the TGA.



Product Number: 300-880-B




Product Number: 300-880-I

Fits bipolar pencils and USA-style two pin forceps - Trewavis Surgical, Valleylab, Conmed, Olsen, Kirwan & Codman.




# STOP

## The Misconnection of Electrosurgical Bipolar Electrodes



**Figure 1. Flexible bipolar electrode cables can allow misconnection of bipolar leads into the monopolar jack, resulting in serious patient injury.**

Source: MDSR: Hazard [Health Devices Jan 1995;24(1):34-5]



"Flying Leads" as shown in Figure 1 and above are now prohibited.

## Trewavis Surgical has the solution!


### ESU Plug End



Trewavis Twin Pin Bipolar Lead

To suit Valleylab or similar twin pin spacing

### Instrument Plug End



Product Number: 300-880-B













Product Number: 300-880-I

Fits bipolar pencils and USA-style two pin forceps - Trewavis Surgical, Valleylab, Conmed, Olsen, Kirwan & Codman.

For more information call Trewavis Surgical or refer to TGA News Issue 56 (July 2008)

# BIPOLAR LEADS

ESU Plug End	Instrument Plug End
 <p>Valleylab Twin Pin <span style="float: right;">300-880-B</span></p>	 <p>Fits bipolar pencils and USA-style two pin forceps - Trewavis Surgical, Valleylab, Conmed, Olsen, Kirwan, Codman.</p>
 <p>Valleylab Twin Pin <span style="float: right;">300-880-I</span></p>	 <p>Same as 300-880-B but with 45° angled plug</p>
 <p>Valleylab Twin Pin <span style="float: right;">300-880-J Heavy Duty Silicon Cable</span></p>	 <p>Fits bipolar pencils and USA-style two pin forceps - Trewavis Surgical, Valleylab, Conmed, Olsen, Kirwan, Codman. With Safety Hood 5 Metre Version available (300-880-N)</p>
 <p>Valleylab Twin Pin <span style="float: right;">300-880-K</span></p>	 <p>Fits Boston Scientific Gold Probe</p>
 <p>Valleylab Twin Pin <span style="float: right;">350-870-C</span></p>	 <p>Fits European flat pin forceps</p>