


DOUBLE INSTRUMENT HOLDERS

For Ring-Handled Instruments



- Designed to meet the new requirements of AS/NZS 4187:2003
- Designed to suit Autoplas sterilization trays (will also fit other commonly used trays)
- 3 Sizes available – 23.0cm Part No. 500-002-D
15.0cm Part No. 500-002-J
11.5cm Part No. 500-002-F
- At a glance, you can see that instruments are not on 1st ratchet
- Easy cleaning
- Dishwasher proof
- Instruments do not tangle, and large sets are easily handled by one person

Spare Clips

Part No. 500-002-E	Single Clip to suit 500-002-D and 500-002-J
Part No. 500-002-H	Packet of 10 to suit 500-002-D and 500-002-J
Part No. 500-002-G	Single Clip to suit 500-002-F
Part No. 500-002-I	Packet of 10 to suit 500-002-F