


## ENT ELECTRODES



300-440-A Scale 1:1



300-440-B Scale 1:1



300-440-C Not to Scale

PRODUCT NUMBER	DESCRIPTION
300-440-A	25mm (1") Insulated Blade on 50mm (2") shaft, 2mm exposed tip
300-440-B	25mm (1") Insulated Blade on 50mm (2") shaft, 5mm exposed tip
300-440-D	25mm (1") Insulated Blade on 100mm (4") shaft, 2mm exposed tip
300-440-F	25mm (1") Insulated Blade on 125mm (5") shaft, 4mm exposed tip
300-440-C	Abbey Sub-Mucous Resection Needle, (Adult Version) 3/16" push pull fitting, lead to suit: 300-730-A
300-440-E	Abbey Sub-Mucous Resection Needle (Paediatric Version) 3/16" push pull fitting, lead to suit: 300-730-A

Electrode shank diameter is standard 2.38mm (3/32")

**Electrodes are also available with a non-stick coating.**

**To order these, add an "N" to the end of the Product Number e.g. 300-440-A would become 300-440-AN.**

