

FUKUDA RETRACTOR




Not to Scale


Used for retracting the humerus to expose the entire glenoid surface.

PRODUCT NUMBER	OVERALL LENGTH	BLADE WIDTH	OPENING
500-001-D	215mm	32mm	26 x 42mm
500-001-E	215mm	38mm	31 x 42mm

HASSON "S" RECTRACTOR



Small



Large

Not to Scale

Used to facilitate the insertion of the Hasson obturator/cannula by retraction of the entry site.
Sold in pairs

PRODUCT NUMBER	DESCRIPTION
500-001-L	"S" Retractor Small
500-001-M	"S" Retractor Large

