


KEON COHEN HEAVY DUTY RETRACTORS

(Modified Gelpi Retractor)



500-001-A
(Left Hand)

500-001-B
(Right Hand)

Not to Scale

Recommended for use in almost any situation where the spinal canal is to be explored.
Also suitable for use in shoulder, knee and hip surgery.

Extended Keon Cohen Retractors (500-001-F, 500-001-G, 500-001-J, 500-001-K) suitable for bariatric patients.

PRODUCT NUMBER	DESCRIPTION	LENGTH	OPENING WIDTH	REACH / DEPTH
500-001-A	Left Hand Standard Reach	200mm	up to 90mm	50mm
500-001-B	Right Hand Standard Reach	200mm	up to 90mm	50mm
500-001-F	Left Hand Extended Reach	200mm	up to 90mm	80mm
500-001-G	Right Hand Extended Reach	200mm	up to 90mm	80mm
500-001-H	Left Hand Standard Reach with bone spike	200mm	up to 90mm	50mm
500-001-I	Right Hand Standard Reach with bone spike	200mm	up to 90mm	50mm
500-001-J	Left Hand Extended Reach with bone spike	200mm	up to 90mm	80mm
500-001-K	Right Hand Extended Reach with bone spike	200mm	up to 90mm	80mm