

AUSTRALIAN MADE SPECIAL MONOPOLAR FORCEPS

MONOPOLAR FORCEPS FOR BREAST SUGERY



PRODUCT NUMBER	DESCRIPTION
300-007	MONOPOLAR FORCEPS 200mm length

300-007 Forceps were developed originally by Trewavis Surgical as a custom made special for Breast Surgery, but can be adapted for use in various surgeries. This instrument can perform the actions of a forcep as well as a blade electrode. The forcep tips are Teflon coated and are designed so that when closed they replicate the action of a blade electrode.

The advantages of using the 300-007 in comparison to a switch pencil with blade electrode and monopolar forceps are:

- Ability to be able to cut and seal vessels with the one instrument
- Avoid tangling of cords between multiple instruments
- Time saving
- Allows Surgeon to react faster in comparison to using 2 instruments
- Reduce blood loss of the patient

Lead to Suit: 300-730-A