



























MONOPOLAR LEADS

ESU Plug End	Instrument Plug End	
 <p data-bbox="276 1093 611 1144">Bovie jack (7.94mm / 5/16") fits the following ESUs:</p> <p data-bbox="276 1171 456 1424">Bard Birtcher Bovie Bowa Concept Conmed/Aspenlab Eschmann Medtrex Neomed Valleylab</p>	 300-710-A	<u>Screw Fitting</u> Forceps, chuckhandles, lap probes with female screw fitting.
	 300-720-A	<u>2.5mm Pin Fitting</u> Forceps, chuckhandles, lap probes with 2.5mm female pin fitting.
	 300-720-B	<u>4mm Male Banana Fitting</u> Forceps, lap probes with 4mm female banana fitting.
	 300-720-C	<u>4mm Recessed Male Banana Fitting</u> Forceps, lap probes with 4mm female banana fitting
	 300-720-D	<u>2mm Pin Fitting</u> Wolf lap cholecystectomy hook electrodes.
	 300-730-A	<u>4.76mm (3/16") Push/Pull Fitting</u> Forceps, chuckhandles, lap probes, with 4.76mm (3/16") male push/pull fitting.
	 300-730-B	<u>4mm Push/Pull Fitting</u> Forceps, chuckhandles, lap probes, lap handles with solid 4mm male pin fitting.
	 300-730-E	<u>4mm Female Banana Fitting</u> Lap probes with 4mm male banana fitting.
	 300-730-F	<u>4mm Stepped Push/Pull Fitting</u> Corson needles, Bugbee electrodes.
	 300-730-G	<u>Light Weight 4mm Push/Pull Fitting</u> Lap handles and instruments with 4mm male banana fitting.
	 300-730-H	<u>Light Weight 4mm Push/Pull Fitting</u> Wolf lap cholecystectomy hook electrodes, Wolf lap handles.
 300-730-I	<u>2.8mm Push/Pull Fitting</u> Cook, Microvasive, AcuSnare Snares.	

All leads have 3.1m (10') silicone cables. (Also available with 5m cables – to order, change 6th digit of Product No. from 0 to 5 e.g. 300-710-A would become 300-715-A for 5m cable.)



MONOPOLAR LEADS

ESU Plug End	Instrument Plug End	
 <p>4mm pin fits the following ESUs:</p> <p>Aesculap Berchtold Erbe T.series Eschmann TD300, TD311, DS302 Martin Sutter BM780</p>	 300-710-A4	<u>Screw Fitting</u> Forceps, chuckhandles, lap probes with female screw fitting.
	 300-720-A4	<u>2.5mm Pin Fitting</u> Forceps, chuckhandles, lap probes with 2.5mm female pin fitting.
	 300-720-B4	<u>4mm Male Banana Fitting</u> Forceps, lap probes with 4mm female banana fitting.
	 300-720-C4	<u>4mm Recessed Male Banana Fitting</u> Forceps, lap probes with 4mm female banana fitting
	 300-720-D4	<u>2mm Pin Fitting</u> Wolf lap cholecystectomy hook electrodes.
	 300-730-A4	<u>4.76mm (3/16") Push/Pull Fitting</u> Forceps, chuckhandles, lap probes, with 4.76mm (3/16") male push/pull fitting.
	 300-730-B4	<u>4mm Push/Pull Fitting</u> Forceps, chuckhandles, lap probes, lap handles with solid 4mm male pin fitting.
	 300-730-E4	<u>4mm Female Banana Fitting</u> Lap probes with 4mm male banana fitting.
	 300-730-F4	<u>4mm Stepped Push/Pull Fitting</u> Corson needles, Bugbee electrodes.
	 300-730-G4	<u>Light Weight 4mm Push/Pull Fitting</u> Lap handles and instruments with 4mm male banana fitting.
	 300-730-H4	<u>Light Weight 4mm Push/Pull Fitting</u> Wolf lap cholecystectomy hook electrodes, Wolf lap handles.
	 300-730-I4	<u>2.8mm Push/Pull Fitting</u> Cook, Microvasive, AcuSnare Snares.

All leads have 3.1m (10') silicone cables. (Also available with 5m cables – to order, change 6th digit of Product No. from 0 to 5 e.g. 300-710-A4 would become 300-715-A4 for 5m cable.)

