

## FUKUDA RETRACTOR




Not to Scale

Used for retracting the humerus to expose the entire glenoid surface.


PRODUCT NUMBER	OVERALL LENGTH	BLADE WIDTH	OPENING
500-001-D	215mm	32mm	26 x 42mm
500-001-E	215mm	38mm	31 x 42mm

## HASSON "S" RECTRACTOR



Small

Not to Scale



Large

Used to facilitate the insertion of the Hasson obturator/cannula by retraction of the entry site. Sold in pairs

PRODUCT NUMBER	DESCRIPTION
500-001-L	"S" Retractor Small
500-001-M	"S" Retractor Large

## FITZPATRICK RECTRACTORS



Not to Scale

Used for Mouth Surgery

PRODUCT NUMBER	DESCRIPTION
500-001-0	Set of 4 Retractors