



REUSABLE LLETZ/LEEP ELECTRODES



Dimensional Example



PRODUCT NUMBER	CROSS BAR LENGTH (A)	LOOP HEIGHT (B)	SHAFT LENGTH (C)	WIRE/BALL TYPE & DIAMETER
300-610-B	10mm	10mm	85mm	Stainless Steel (0.3mm)
300-610-C	20mm	15mm	85mm	Stainless Steel (0.3mm)
300-610-D	25mm	15mm	85mm	Stainless Steel (0.3mm)
300-610-E		4mm Ball	85mm	Stainless Steel
300-610-F	15mm	15mm	85mm	Stainless Steel (0.3mm)
300-610-H	25mm	20mm	85mm	Stainless Steel (0.3mm)
300-611-B	10mm	10mm	130mm	Stainless Steel (0.3mm)
300-611-C	20mm	15mm	130mm	Stainless Steel (0.3mm)
300-611-D	25mm	15mm	130mm	Stainless Steel (0.3mm)
300-611-E		4mm Ball	130mm	Stainless Steel
300-612-B	10mm	10mm	85mm	Tungsten (0.2mm)
300-612-C	20mm	15mm	85mm	Tungsten (0.2mm)
300-612-D	25mm	15mm	85mm	Tungsten (0.2mm)
300-613-B	10mm	10mm	130mm	Tungsten (0.2mm)
300-613-C	20mm	15mm	130mm	Tungsten (0.2mm)
300-613-D	25mm	15mm	130mm	Tungsten (0.2mm)
300-613-F	15mm	15mm	130mm	Tungsten (0.2mm)
300-613-G	20mm	20mm	130mm	Tungsten (0.2mm)
300-613-H	25mm	20mm	130mm	Tungsten (0.2mm)
300-615-B	12mm	8mm	100mm	Stainless Steel (0.23mm)
300-615-C	15mm	10mm	100mm	Stainless Steel (0.23mm)
300-615-D	20mm	12mm	100mm	Stainless Steel (0.23mm)
300-616-B	12mm	8mm	130mm	Stainless Steel (0.23mm)
300-616-C	15mm	10mm	130mm	Stainless Steel (0.23mm)
300-616-D	20mm	12mm	130mm	Stainless Steel (0.23mm)
300-617-B	12mm	8mm	100mm	Tungsten (0.2mm)
300-617-C	15mm	10mm	100mm	Tungsten (0.2mm)
300-617-D	20mm	12mm	100mm	Tungsten (0.2mm)
300-618-B	12mm	8mm	130mm	Tungsten (0.2mm)
300-618-C	15mm	10mm	130mm	Tungsten (0.2mm)
300-618-D	20mm	12mm	130mm	Tungsten (0.2mm)

Electrode shank diameter is standard 2.38mm (3/32"), Ball Electrodes available with a non-stick coating.