

SINGLE USE STERILE LLETZ/LEEP ELECTRODES


LOOP ELECTRODES

PRODUCT NUMBER	CROSS BAR LENGTH	LOOP HEIGHT
TSS001	10 mm	10 mm
TSS002	12 mm	8 mm
TSS003	15 mm	10 mm
TSS004	15 mm	15 mm
TSS005	20 mm	8 mm
TSS006	20 mm	12 mm
TSS007	20 mm	15 mm
TSS008	25 mm	20 mm




BALL ELECTRODES

PRODUCT NUMBER	BALL SIZE	
TSS011	3 mm	
TSS012	5 mm	



NEEDLE ELECTRODES

PRODUCT NUMBER	NEEDLE LENGTH	EXPOSED TIP
TSS013	16 mm	3 mm



- 120mm long shaft, standard 2.38mm (3/32") diameter
- High grade tungsten wire Loop Electrodes
- Stainless Steel Ball and Needle Electrodes
- Packed 5 per pack

