

# TUNGSTEN MICRO DISSECTION NEEDLE ELECTRODES



300-470-A Scale 1:1

PRODUCT NUMBER	DESCRIPTION
300-470-A	30mm Straight Needle, 3mm exposed tip
300-470-B	20mm Straight Needle, 3mm exposed tip
300-470-C	40mm Straight Needle, 3mm exposed tip
300-470-D	30mm Needle, 3mm exposed tip at 45 degree angle
300-470-E	112mm Straight Needle, 3mm exposed tip

Electrode shank diameter is standard 2.38mm (3/32")

# STERILE SINGLE USE TUNGSTEN MICRO DISSECTION NEEDLE ELECTRODES



PRODUCT NUMBER	DESCRIPTION
TSS020	50mm Straight Needle, Single Use, Sterile
TSS021	60mm Straight Needle, Single Use, Sterile

- Tungsten Needle causes Minimal Tissue Injury
- Reduced Scarring and Better Wound Healing
- Packed as 5 individual sterile units per Pack
- Sterile, Ready for Use

